

MULTIGEN⁺

THE TOTAL KNEE SYSTEM

EVOLUTE INSTRUMENT SET

Evolute Instrument Set:

- Precision design
- Compact instrumentation layout
- Intuitive technique from primary through to revision
- Smaller soft tissue friendly blocks



EVOLUTE Instrument Set

The EVOLUTE Instrument Set is designed to ensure a reproducible surgical technique which is clear and simple. This ensures a well positioned implant for every patient, whilst allowing for minimal disruption to the soft tissue mechanism.

These ingenious instruments are straightforward for you and your OR team, ensuring a precise, intuitive technique time and time again. They incorporate insightful and unique concepts, whilst allowing a versatile approach to knee replacement, paving the way for every surgeon's preferred surgical approach. The same basic instrument set is utilised for the entire **MULTIGEN⁺ Total Knee System**.

Unique features have been incorporated into the knee system, designed by surgeons and engineers who understand the challenges of knee surgery. They have taken these complexities and made them easier to deal with, through one basic instrument set.

With precision design, compact instrumentation layout and soft tissue friendly cutting blocks, the EVOLUTE Instrument Set offers intuitive primary surgery.

The Femur

Modular Distal Femoral Mask Guide:

- Three-part modular instrument for quick and easy assembly, even within the joint space
- Low profile ensures minimum soft tissue disruption
- Precise distal femoral cut to allow for a well balanced knee in extension

Femoral Calliper Device:

- Accurate sizing prior to any bone cuts
- Magnetic modules ensure quick connection and disconnection during surgery
- External rotation provided through 0°-3° and 5° module options with one main calliper device

Four-in-One:

- Low profile for soft tissue friendly surgery
- Stable fixation on the femur during bone cuts
- Clearly marked sizes for reference
- Optimized blade capturing slots to ensure accurate bone cuts



The Tibia

Low-Profile Tibial Blocks:

- MIS block geometry allows clear visibility during resection
- Designed to protect the soft tissues
- Standard block geometry
- Fully adjustable without the need to re-pin

Proximal Resection Guide:

- Fully adjustable to ensure accurate initial positioning
- Allows medial and lateral reference of the tibial plateau
- Allows adjustable posterior slope
- Three point fastening ensures stability throughout the resection
- Secure locking wheels

Instrument Configuration

The EVOLUTE Instrument Set includes a unique compact tray layout and additional features to reduce the number of instruments, whilst providing an intuitive technique for you and your surgical team.

B.6610.83.020.1

081000

